

Eye protection/Eyeshields

Note: The information given in this catalogue regarding personal safety equipment is of general nature and does not indicate whether the product is appropriate to the safety-requirements of a specific workplace. Please check the information given by the manufacturer in detail in case of any doubts.

European standards for eye protection

Basic standards

EN 166	Guidelines
EN 167	Test methods for optical glasses
EN 168	Test methods for non-optical procedures

Standards acc. to usage category

EN 170	UV filter
EN 171	IR filter
EN 172	UV protection for industrial applications

1



1 LLG-Safety eyeshields, clear

NEW!

- CE approved to EN166 and EN170
- Visitor Spec with single polycarbonate lens
- Can be worn over prescription glasses
- Side venting
- Lense: clear
- Weight: 43g

Colour	Lens	PK	Cat. No.
clear	clear, uncoated	10	9.006 002

2



2 LLG-Safety eyeshields, blue

NEW!

- CE approved to EN166 and EN 170
- Ultra-light allergic tested frame
- Integrated sideshields
- Temple adjustable in length and inclination
- Weight: 34g
- Frame colour: blue
- Lense: clear
- Anti-scratch lense

Colour	Lens	PK	Cat. No.
blue	clear, anti-scratch lens	1	9.006 001

3



3 LLG-Safety eyeshields, black/orange

NEW!

- CE approved to EN166 and EN170
- Model with sporty and colourful line
- High wearability
- Over moulded nose bridge and temples for maximum comfort
- allergenic rubber eyebrow protection
- Temple adjustable in length
- Weight: 30g
- Frame colour: black/orange
- Lense: clear
- Anti-scratch and anti-fog
- Protection against UV400

Colour	Lens	PK	Cat. No.
black/orange	clear, anti-scratch and anti-fog lens, UV protection	1	9.006 003

4



9.005 072

Safety eyeshields CLAREX

- DIN EN 166 1 - FT
- sports-style design
- adapts to any head shape
- excellent all-round view
- scratch-resistant and anti-fog lens
- 99.99% UV protection using innovative UV 380 technology
- 4-position temple length adjustment

EKASTU Safety

5



9.005 071

Type	PK	Cat. No.
CLAREX, transparent	1	9.005 072 4
CLAREX, tinted	1	9.005 071 5

1 Safety eyeshields x-trend 9177

With an anatomically shaped, wrap-around lens, the uvex x-trend is an expression of today's fashion - after all, eyewear with large lenses conquered the world long ago. In addition to its stylish design, the uvex x-trend offers a high degree of safety with its large lens. It offers a clear view in all directions, 20 % more cover than traditional spectacles, gives UV protection of up to 400 nm and comfortable temples.

Coating of the lens ultradura: excellent scratch resistance, 100 % UV protection

Colour	Lens	PK	Cat. No.
black	PC transparent / UV 2-1.2, ultradura	1	9.005 147
black	PC grey / UV 5-2.5, ultradura	1	9.005 149

1

Uvex



2 Safety eyeshields Skylite 9174

Red fashionable sports-style frame. Arms can be tilted. High impact resistant, optically correct, polycarbonate lens. Lenses have integral matt side panels. Weighing only 26g, these eyeshields are the lightest of their kind.

The lens coating is ultradura, excellent scratch resistance and gives 100% UV protection.

Colour	Lens	PK	Cat. No.
red	grey / UV 5-2.5 ultradura	1	9.005 127
red	transparent / UV 2-1.2 ultradura	1	9.005 128
blue	grey / UV 5-2.5 ultradura	1	6.200 303
blue	transparent / UV 2-1.2 ultradura	1	6.900 098

2

Uvex



3 Safety eyeshields super g 9172

The lightest safety eyeshields worldwide.

Using state-of-the-art manufacturing processes, an extremely thin lens with exceptional optical clarity can be created. The lenses have ideal aerodynamics and are angled forwards to provide excellent ventilation, keeping the eyes comfortable.

The frame is made of ultra light, high-tech polymers, guaranteeing maximum flexibility. The innovative side piece concept features Uvex x-stream-technology for an excellent, non-slip fit and a wide range of applications.

The super g is completely metal-free and provides top-level protection, with a weight of only 18g.

The coating of the lens is supravision HC-AF, outside scratch-resistant, inside non-fogging and gives 100% UV protection. Optidur NC, excellent scratch and chemical resistance and gives 100% UV protection.

Colour	Lens	PK	Cat. No.
crystal	PC transparent/UV, 2-1.2, supravision HC-AF	1	6.228 096
titan	PC transparent/UV, 2-1.2, optidur NC	1	9.005 164
navy blue	PC transparent/UV, 2-1.2, optidur NC	1	9.005 169

3

Uvex



4 Eyeshields i-vo 9160

Fashionable, sports-style design, eyeshield with comfortable, semi-elastic frame incorporating two-component technology. Ultra-transparent lenses with scratch-proof, outer coating. Eyeshields for use when carrying out all potentially hazardous work in the laboratory.

The coating of the lens is supravision HC-AF, outside scratch-resistant, inside non-fogging and gives 100% UV protection. Optidur NC, easy to clean coating, excellent scratch and chemical resistance and gives 100% UV protection.

Colour	Lens	PK	Cat. No.
blue/grey	transparent/UV 2-1.2, supravision HC-AF	1	9.005 116
blue/orange	transparent/UV 2-1.2, supravision HC-AF	1	6.311 677
black/grey	transparent/UV 2-1.2, supravision HC-AF	1	6.901 653
blue/orange	transparent/UV 2-1.2, optidur NC	1	6.324 668
black/grey	grey/UV 5-2.5, optidur NC	1	9.005 117
blue/orange	brown/UV 5-2.5, optidur NC	1	6.318 848

NEW!

4

Uvex



Eye protection/Eyeshields

1


1 Spectacles uvex racer MT 9153

- colour: titan
- metal sunglasses, sports-type
- very low-profile metal frame, high flexibility and extremely low weight
- lens held in place without screws, making the glasses safe, even under mechanical strain, the glasses fit perfectly on all head sizes
- flexible nose piece, adjustable
- scratch-resistant lens

Uvex

Colour	Lens	PK	Cat. No.
gunmetal grey	transparent/UV 2-1.2, non-fogging	1	9.005 134

2


2 Industrial safety eyeshield Pulsafe Millennia

Pulsafe Millennia offers uninterrupted vision through 9 dioptre lens curves with constant distance from pupil, close fitting to the eye sockets, ergonomically shaped, non-pressing frame in a choice of colours, and modern styling. With non-misting FogBan lenses. Even in extreme situations such as erratic changes in temperature, high humidity etc. condensation is no longer an issue.

Colour	Lens	PK	Cat. No.
Black	Transparent, non-misting	1	9.005 073
Black	Grey, non-misting	1	9.005 074
Blue	Transparent, non-misting	1	9.005 075

3


3 Superlight safety eyeshields Uvex Astrospec 9168

- Unrestricted peripheral vision
- Weight 36g
- With Uvex Duo-flex® frame for maximum comfort
- 3-position adjustable lens
- 4 adjustable side-arm lengths with built-in pull-out stops
- Lens made from impact resistant polycarbonate, 100% UV protection
- Integral side shields

Uvex

Coating of the lens
 optidur NC: excellent scratch and chemical resistance
 optidur 4C PLUS: anti-fog lens and scratch resistant

Colour	Lens	PK	Cat. No.
Silver-grey	transparent / 2-1.2 UV, optidur NC	1	9.005 105
Yellow/black	transparent / 2-1.2 UV, optidur NC	1	9.005 106
Blue	transparent / 2-1.2 UV, optidur 4C PLUS	1	9.005 109
Yellow/black	grey / 5-2.5 UV, optidur NC	1	9.005 100
Violet/green	transparent / 2-1.2 UV, optidur NC	1	9.005 157

4


4 Eyeshield astroflex 9163

Highly suited to overhead work. The soft, flexible sealing rim protects the face from frame pressure, dust and splashes from above. Additionally provides impact protection.

Uvex

Coating of the lens supravision HC-AF:
 outside scratch resistant, inside non-fogging, 100% UV protection

Colour	Lens	PK	Cat. No.
blue	transparent/UV 2-1.2, supravision HC-AF	1	9.005 118



1 Safety eyeshields Skyper 9195/Skyper S 9196

colour: blue, skyper 9195

Uvex

Extremely lightweight (approx. 30g), modern, safety eyeshields with panoramic field of view and excellent fit. Five individual side-arm lengths with integral pull-out stops. The soft frame ensures a pressure-free fit and can be finely adjusted. Lenses and frame are made of impact resistant, synthetic material. Matt finish side shields offer extra protection from reflected light

colour: blue, skyper S 9196

Skyper S is a smaller model of Skyper 9195, specially designed for smaller head sizes. Provide the same features as Skyper 9195.

Coating of the lens

supravision HC-AF: outside scratch-resistant, inside non-fogging, 100 % UV protection

optidur NC: excellent scratch and chemical resistance, 100 % UV protection

Type	Lens	PK	Cat. No.
Skyper 9195	transparent/UV, 2-1.2, supravision HC-AF	1	9.005 123
Skyper 9195	transparent/UV, 2-1.2, optidur NC	1	9.005 125
Skyper S 9196	transparent/UV, 2-1.2, supravision HC-AF	1	9.005 124
Skyper S 9196	transparent/UV, 2-1.2, optidur NC	1	9.005 126



2 Safety eyeshields JACKSON SAFETY*V60 Nemesis RX

NEW!

Kimberly-Clark

A stylish new standard in eye protection that improves personal comfort and safety and assists worker performance.

- Impact resistant polycarbonate lens that meets European standard EN166 1F
- Lightweight wrap-around design for added comfort and coverage
- 99.9% UVA/UVB Protection
- Black frame
- Bifocal style with diopters for vision assistance
- V60 Nemesis RX is a stylish solution for users who require protection reading diopters (strengths include), +1.0 to +3.0 available
- Enhanced nose piece design to channel away sweat
- Includes free neck cord

Diopter	PK	Cat. No.
+1.0	1	9.005 171
+1.5	1	9.005 172
+2.0	1	9.005 173
+2.5	1	9.005 174
+3.0	1	9.005 175



3 Safety eyeshields skyguard NT 9175

NEW!

Uvex

The uvex skyguard NT spectacles have been improved to provide all-round protection.

Based on the existing uvex skylite spectacle we have added a conceptually simple, injection moulded TPU frame to significantly increase wearer protection levels, through minimising the risk of particle ingress. The flexible, closed-pore soft components of the frame enables a customised, individual fit and are easy to clean to promote long term use. Optional: uvex side arms can be adjusted to four different lengths and finely angled, guaranteeing a comfortable, customised fit.

Colour	Lens	PK	Cat. No.
blue/grey	clear/UV 2-1.2 supravision extrem	1	9.005 078
grey/orange	clear/UV 2-1.2 supravision HC-AF	1	6.239 420



4 Panoramic eyeshield ultrasonic 9302

Extremely safe and comfortable, semi-elastic frame incorporating two-component technology. Panoramic lens with two separate lacquer coatings.

Uvex

Coating of the lens

supravision HC-AF: outside scratch-resistant, inside non-fogging, 100 % UV protection

Colour	Lens	PK	Cat. No.
orange/grey	transparent / UV 2-1.2, supravision HC-AF	1	9.005 115



Spare lens for panoramic eyeshield ultrasonic 9302

Coating of the lens supravision HC-AF: outside scratch-resistant, inside non-fogging.

Uvex

Description	PK	Cat. No.
Spare lens for panoramic eyeshield ultrasonic 9302	1	6.401 441

1



1 Panoramic eyeshield carbonvision 9307

NEW!

Uvex

The intelligent use of ultra modern materials and manufacturing technologies enables uvex to design frames which combine reduced size/volume and high mechanical resistance. uvex carbonvision provides the wearer with unparalleled comfort levels despite the very high level of mechanical stability. The TPU layer is permanently fixed to the frame, ensuring a perfect individual fit and its easy to clean and maintain. The uvex carbonvision, weighing just 46g is undoubtedly one of the lightest wide-vision goggles on the market, enabling the wearer to work for longer periods without the distractions and problems caused by discomfort, thereby increasing his or her productivity. The newly developed supravision extreme coating technology prevent the goggles from steaming up, even in the most demanding conditions.

Uvex carbonvision goggles are recommended for protection against dust, liquid splashes and mechanical hazards. The goggles are ideal for wearing with disposable respirators and bump caps.

Colour	Lens	PK	Cat. No.
black/grey	clear/UV 2-1.2 uvex supravision extreme	1	9.005 079

2



2 Panoramic eyeshield uvex 9405

Uvex

- large goggle with front ventilation system
- suitable for use with half masks and filtering half masks
- ideal for people who wear prescription glasses, easy headband adjustment, easy changing of lens
- crystal acrylic lens, non-fogging

Colour	Lens	PK	Cat. No.
transparent grey	transparent /UV 2-1.2, non-fogging	1	9.005 090

3



3 LLG-Visitor safety eyeshields, clear

NEW!

- CE approved to EN166 and EN170
- Visitor Spec with single polycarbonate lens
- Can be worn over prescription glasses
- Side venting
- Lense: clear
- Weight: 43g

Colour	Lens	PK	Cat. No.
clear	clear, uncoated	10	9.006 002

4



4 Eyeshield Armamax AX1H

- Made of tough-coated PC.
- Covers the whole eye area
 - Complete side protection with wide frame
 - Matt finish browpiece prevents glare
 - High nosebridge prevents scratches on the lens
 - Extra large
 - Fits over almost all prescription spectacles

Type	PK	Cat. No.
Eyeshield	1	9.005 065

5



5 Overgoggles Type 9161 and 9161 duo-flex®

Uvex

- designed to fit over prescription spectacles
- unlimited side vision through super panorama lenses with a 180° field of view
- excellent, pressure-free fit
- effective indirect ventilation through special frame construction with side shields

Type 9161 duo-flex®

- frame can be adjusted four ways for optimal comfort
- scratch resistant lens

Type	Colour	Lens	PK	Cat. No.
9161	transparent*	transparent, PC lenses	1	9.005 166
9161 duo-flex®	blue/black	transparent, UV 2-1.2 optidur 3000	1	9.005 107

*autoclavable

1 Panoramic vision safety goggles super OTG 9169

NEW!

Uvex

- does not pinch, fits perfectly
- lightweight eyeshield with unrestricted peripheral vision
- protective shielding thanks to flat base curve lens
- lens inclination of 11° ensures good downwards view
- x-stream-technology (XST), flexible temple concept for exceptional user-comfort, with or without prescription glasses
- highly scratch-proof lens

Coating of the lens optidur NC:
excellent scratch and chemical resistance, 100% UV protection

Colour	Lens	PK	Cat. No.
navy blue	PC grey/UV 2-1.2, optidur NC	1	6.234 320
navy blue	PC transparent/UV 2-1.2, optidur NC	1	7.650 112



2 Safety spectacles dispenser

Guards against dust.
Hinged cover.
With holes for wall mounting and label in German.
To hold single pair of spectacles.

Internal dimensions (W x D x H) mm	PK	Cat. No.
180 x 85 x 75	1	9.005 131



3 Safety Glasses Dispenser Clearly Safe®

NEW!

Heathrow Scientific

Ergonomically angled for easy dispensing. Translucent acrylic dispenser provides easy access to safety accessories. Flip open to fill from top; close to keep contents dust-free. Tilt-open bottom lid stays in place during repeated use. The safety glasses dispenser holds approximately 20 glasses.

Description	Dimensions (W x D x H) mm	PK	Cat. No.
Clearly Safe®, Safety Glasses Dispenser	203 x 203 x 406	1	6.238 548



4 5 SecuBox Midi

Necessary for safe and sanitary storage of protective equipment at work.

Overall dimensions: 23.6 x 22.5 x 12.5 cm
Internal dimensions: 23.0 x 21.8 x 11.8 cm

Colour	Type	PK	Cat. No.
blue	Eye protection, without inlet	1	9.005 500
transparent	Eye protection, without inlet	1	9.005 501
blue	Ear protection	1	9.005 502
transparent	Ear protection	1	9.005 503
blue	Eye + ear protection	1	9.005 504
transparent	Eye + ear protection	1	9.005 505
blue	Breathing protection	1	9.005 506
transparent	Breathing protection	1	9.005 507
blue	Protection gloves	1	9.005 508
transparent	Protection gloves	1	9.005 509
blue	Protection gloves, long	1	9.005 511
transparent	Protection gloves, long	1	9.005 512
blue	Hair protection	1	9.005 513
transparent	Hair protection	1	9.005 514
blue	Protective overall	1	9.005 515
transparent	Protective overall	1	9.005 516
anthrazit	Inlet for Box	1	9.005 531



1 2 SecuBox Maxi

Necessary for safe and sanitary storage of protective equipment at work.

Overall dimensions: 23.6 x 31.5 x 20.0cm

Internal dimensions: 23.0 x 29.9 x 19.2cm

Colour	Type	PK	Cat. No.
blue	Eye protection with insert for 12 safety glasses	1	9.005 549
blue	Eye and ear protection	1	9.005 544
blue	Face shield	1	9.005 547
blue	Helmet	1	9.005 541
blue	Helmet and eye protection	1	9.005 542
blue	Helmet and ear protection	1	9.005 543
blue	Helmet and full face mask	1	9.005 546
blue	Full face mask	1	9.005 545



3 Cleaning tissues 9963

- suitable for all uvex lenses
 - silicone-free
 - 100 individually sealed tissues in a box.
 - incl. clip for wall mounting
- Please order wall mounting unit separately.

Uvex

Description	PK	Cat. No.
Cleaning tissues 9963	100	6.200 654



4 Lens cleaning station 99700

For wall mounting. Easy operation. Frequent cleaning increases safety and enhances durability of spectacles.
Refills: Silicone-free paper 9.005 102, cleaning fluid 9.005 103 and plastic pump 9.005 104

Uvex

Description	PK	Cat. No.
Station, complete, with paper, plastic pump bottle and cleaning fluid	1	9.005 101

Accessories for lens cleaning station 99700

Uvex

Description	PK	Cat. No.
Silicone-free paper	1	9.005 102
Cleaning fluid, 0.5 L	1	9.005 103
Plastic pump	1	9.005 104

1 Face shields Clearways



1



Browguard CB14:

- light weight
- elastic napestrap allows individual adjustment
- adjust headband while wearing the face shield

Browguard CB20:

- light weight
- browguard with ratched headband

Face shield visors:

- impact and chemical resistant
- comfortable and lightweight, even during extended usage
- visor can be tilted upwards and locked in any position
- easy exchangeable visors

Please order browguard and visor separately.

Type	Description	PK	Cat. No.
CB14	Browguard with elastic napestrap	1	9.005 067
CB20	Browguard with ratched headband	1	7.617 935
CV83 P	Face shield visor, PC, against mechanical hazards	1	9.005 066
CV84 A	Face shield visor, Acetat, against spilling chemical	1	6.070 565

2 Protective face shield CONTRACID I



2



Shield headband made from vulcanised fibre. With interchangeable cellulose acetate visor. The large visor (410 x 285mm) gives unrestricted vision and also protects areas below and above the face. Visor can be tilted up and locked in any position.

Type	PK	Cat. No.
Protective face shield 0.8mm	1	9.005 145
Protective face shield 1.15mm, UV-protection 3-1.2 (EN 170)	1	9.005 146
Replacement panels 0.8mm (PET)	1	9.005 142
Headband for contracid face shield CONTRACID I	1	6.200 935

3 Cryo-temperature protection kit CRYO-TEMP-SHIELD®

Consisting of:

- 1 pair Cryo Gloves®, size: optional, **see table**
- 1 x face protection shield with adjustable head mounting
- 2 x safety labels, 300mm diameter, entitled "Face protection" and "Wear gloves"
- 1 x protection apron, Cryo Apron®, width: 610mm, length: 1070mm

Type	Type	Size	PK	Cat. No.
1001-S	Gloves forearm length	S (8)	1	9.405 140
1001-M	Gloves forearm length	M (9)	1	9.405 141
1001-L	Gloves forearm length	L (10)	1	9.405 142
1001-XL	Gloves forearm length	XL (11)	1	9.405 143
1102-M	Gloves elbow length	M (9)	1	9.405 144
1102-L	Gloves elbow length	L (10)	1	9.405 145
1102-XL	Gloves elbow length	XL (11)	1	9.405 146
1103-M	Gloves shoulder length	M (9)	1	9.405 147
1103-L	Gloves shoulder length	L (10)	1	9.405 148
1103-XL	Gloves shoulder length	XL (11)	1	9.405 149

3



Eye protection/Eye wash bottles



Eye-wash bottles, Barikos

Barikos KS

Barikos KS and Mini Barikos KS are maintenance-free models of eye wash bottles. Containing water, they can be kept ready for use for at least two years in the sealed state. Both models are ideal for including within first aid boxes and sets. All above mentioned Barikos eyewash bottles comply with DIN EN 15154-4. A protective cover clips onto the spray head of the eye wash bottle thus preventing any leakage. The tamper-evident seal guarantees proper function.

Barikos K

Barikos and Barikos K are designed to be filled with clear water. If the user is likely to be operating an eyewash bottle when standing, we recommend the Barikos. The user need only incline his head slightly downwards, ensuring good contact with the spray head of the bottle. The Barikos K-eyewash bottle (vol.620 ml) can be used immediately, even when the injured person is lying down. As soon as the bottle is tilted downwards, a patented sliding valve opens and the washing fluid reaches the spray head unhindered. The Mini Barikos eyewash bottle (vol.175ml) is designed for persons who are frequently present in different working areas.

Both bottles are suitable for use with upright and with the downward spray head.

Type	PK	Cat. No.	
Barikos KS, 620ml, filled	1	9.005 009	1
Mini Barikos KS, 175ml, filled	1	9.005 008	2
Barikos, 620ml, empty	1	9.005 011	
Mini Barikos, 175ml, empty	1	9.005 020	
Barikos K, 620ml, empty	1	9.005 014	
Barikos K / KS wall unit for 1 eye wash bottle	1	9.005 015	3
Barikos K / KS wall unit for 2 eye wash bottles	1	9.005 017	4
Plastic cabinet for wall assembly for 1 eye wash bottle	1	6.076 057	5



oculavNIT® Eye Wash

oculav NIT® Eye Wash

First aid eye wash sterile solution. 250ml of sterile solution for eye wash treatment at the accident scene. Especially effective against eye injuries stemming from acid and caustic burns as well as pepper, CS and tear gases.

- easy to dispense, providing gentle removal of foreign bodies
- temperature-resistant, avoids contamination, maintenance-free
- one handed use
- according to MPG (German Act on Medical Devices) DIN EN 15154-4.

oculav NIT®-Box

Orange plastic box, size 260 x 160 x 80mm, with printed inscriptions and practical wall brackets. With 4 ready to use spray bottles oculav NIT® Eye Wash. The first-aid set for eye injuries contains DuOcul eye compresses, eye flaps, plaster strips, compresses and disposable gloves.



Description	PK	Cat. No.	
oculav NIT® Eye Wash, 250ml bottle	1	9.005 001	6
oculav NIT®-Box with 4 x 250ml bottles	1	9.005 002	7

1 | 2 Eye-wash bottle, 0.9% NaCl

- choice of 0.2 litre or 0.5 litre bottles
- eye cup - keeps the eyes open during application
- saline solution (0.9%)
- reuse isn't possible due to sealing
- according to DIN EN 15154-4
- CE-marked
- shelf life: 3 years

Breuell Ingenieurbüro

Applications:

- mobile application
- transportation of dangerous goods acc. GGVS
- emergency medication
- First-aid boxes

Type	PK	Cat. No.
Eye-wash bottle 200 ml	1	9.733 794
Eye-wash bottle 500 ml	1	9.733 793

1



2



3 Eye wash bottle pH-neutral

A small, handy bottle containing 4.9% sterile phosphate buffer solution, that neutralises the eye's fluid after accidents involving acids and alkalies. The bottle is easy to take with you in the specially designed belt bag, the first aid box, the toolbox etc. A wall holder is also available. The bottle is provided with an ergonomic eye cup, a dust cap and clear instructions. For refilling Eyewash station.

Breuell Ingenieurbüro

- durability: 3 years
- flush time: about 2 minutes
- according to DIN EN 15154-4

Description	PK	Cat. No.
Spare eyewash bottle, filled, 200 ml	1	9.733 788

3



4 Eyewash emergency station DUO

NEW!

For parallel treatment of both eyes. With one pH neutral DUO eye wash bottle with 500ml sterile phosphate buffer solution (4.9%) and one DUO eye wash bottle with 1000ml sterile sodium chloride solution (0.9%). With wall mounting, enclosed pictogram and mirror. Inclusive eye emergency plan. Especially suitable for workplaces where chemicals and contaminants can penetrate the eye. Advice: Also after treatment with pH neutral solution continued flushing until medical treatment with sodium chloride solution is recommended. Shelf life: 3 years. According to DIN EN 15154-4:2009.

Breuell Ingenieurbüro

Type	PK	Cat. No.
Eyewash emergency station DUO	1	9.733 816
Eye wash-bottle pH neutral DUO	1	9.733 817
Eye cleaning station DUO incl. 1 bottle NaCl	1	9.733 818

4



5 Eye wash emergency station, wall-mounting

A dust-proof wall mounted station containing one bottle 200ml sterile phosphate buffer solution and one bottle 500ml sterile sodium chloride solution complete with pictogram and a mirror. Especially suitable for workplaces where on the spot rinsing is required, smaller workplaces, and mobile workplaces. According to DIN EN 15154-4.

Breuell Ingenieurbüro

Dimensions of the station (H x W x D): 265 x 227 x 110mm
 Durability: 3 years
 Flush time: about 7 minutes

Type	PK	Cat. No.
Eyewash emergency station, wall-mounting	1	9.733 787
Spare eyewash bottle, pH-neutral, filled, 200ml	1	9.733 788
Spare eyewash bottle, filled, 500ml	1	9.733 793

5



Eye protection/Eye wash stations



9.733 779

3 Hand-held drench hoses

Two sprays in one. The hose can be used as a fixed eyewash when fitted into the mount. This allows the injured person to use both hands to hold open the eyes. The hose can also be removed from the mount so that other body parts and reclining patients can be treated.

Breuell Ingenieurbüro

- For wall and table mounting
- This drench hose meets EN 15154-2:2006 requirements
- DIN certification label included
- The spray heads ensure a wide distribution of water so that areas around the eyes are also covered and the exact positioning of the eyes is therefore not critical
- Large holes in the spray heads protect against limescale. The spray heads can be removed for cleaning purposes
- Rubber shrouds on the nozzles help prevent accidental injury
- Self-closing valve with squeeze lever operation
- Lightweight dust caps fall off automatically on operation
- Integral jet regulator keeps the water jet constant independent of the water pressure
- Built-in non-return valve prevents any backflow of contaminated water



9.733 786

Spray height: 250mm
 Hose length: 1500mm
 Connection: 1/2" swivel nut
 Water output: 7 litres/min per spray head, regulated
 Inlet fitting: M 28 x 80mm.

Type	PK	Cat. No.
1 spray head, 45° angle	1	9.733 778
2 spray heads, 45° angle	1	9.733 779 1
1 spray head, vertical	1	9.733 786 2



4 Wall/table bracket for hand-held eye showers

- table mounting dimensions: M28 x 80mm
- installation dimension: 30mm
- fastening-parts kit: for surface mounting
- seal to table surface: with joint ring

Breuell Ingenieurbüro

Type	PK	Cat. No.
Wall-/table bracket for hand-held eye showers	1	9.733 768



5 Body showers



Special showering heads ensure that large water quantities in gentle streams reaches the injured person without damaged the injured skin.

- The showers meet EN 15154-1:2006
- Marked sign according to DIN is included
- Large dimensioned maintenance-free ball tap for large throughput quantities
- Throughput quantity: approx. 60 litres/min.
- Simple assembly
- Stable structure
- Chemical-resistant PE powder painting
- High speed showering heads made from plastic, with improved spray pattern. Corrosion resistant, calcify- and maintenance free, very durable, even emptying.
- Large, easy-to-use release bar
- Comes with all mounting materials
- Connection 3/4"

Inlet	Attachment	Length mm	PK	Cat. No.
Surface mount	Wall flange	625	1	9.733 810
Concealed mount	Wall/60 mm - 3/4" pipe	625	1	9.733 811
Surface mount	Over door wall flange	565	1	9.733 812