

**1 2 3 4** LLG-GHS warning labels, self-adhesive, roll in dispenser box



Dispenser box of 250 labels. These labels are produced using an extremely resistant 100µ PE film.  
 Due to their special surface coating they are resistant to most solvents and excellent for indoor or outdoor use.  
 The strong adhesive, although permanent, can be removed from glass containers without leaving any trace.  
 Trilingual signal word: German, English, French.  
 Available in 2 sizes.



Type	Description	Label	Dimensions mm	PK	Cat. No.
GHS 01	warning	exploding bomb	26 x 37	250	9.105 700
GHS 01	danger	exploding bomb	26 x 37	250	9.105 701
GHS 02	warning	flame	26 x 37	250	9.105 702
GHS 02	danger	flame	26 x 37	250	9.105 703
GHS 03	warning	flame over circle	26 x 37	250	9.105 704
GHS 03	danger	flame over circle	26 x 37	250	9.105 705
GHS 04	warning	gas cylinder	26 x 37	250	9.105 706
GHS 05	warning	corrosive	26 x 37	250	9.105 707
GHS 05	danger	corrosive	26 x 37	250	9.105 708
GHS 06	danger	skull and crossbones	26 x 37	250	9.105 709
GHS 07	warning	exclamation mark	26 x 37	250	9.105 710
GHS 08	warning	health hazard	26 x 37	250	9.105 711
GHS 08	danger	health hazard	26 x 37	250	9.105 712
GHS 09	warning	environment	26 x 37	250	9.105 713
GHS 01	warning	exploding bomb	37 x 52	250	9.105 720
GHS 01	danger	exploding bomb	37 x 52	250	9.105 721
GHS 02	warning	flame	37 x 52	250	9.105 722
GHS 02	danger	flame	37 x 52	250	9.105 723
GHS 03	warning	flame over circle	37 x 52	250	9.105 724
GHS 03	danger	flame over circle	37 x 52	250	9.105 725
GHS 04	warning	gas cylinder	37 x 52	250	9.105 726
GHS 05	warning	corrosive	37 x 52	250	9.105 727
GHS 05	danger	corrosive	37 x 52	250	9.105 728
GHS 06	danger	skull and crossbones	37 x 52	250	9.105 729
GHS 07	warning	exclamation mark	37 x 52	250	9.105 730
GHS 08	warning	health hazard	37 x 52	250	9.105 731
GHS 08	danger	health hazard	37 x 52	250	9.105 732
GHS 09	warning	environment	37 x 52	250	9.105 733



**5** Hazard and prohibitive labels

Self-adhesive vinyl.

Type	Dia. mm	PK	Cat. No.
Warning: Flammable substances	20	20	9.105 340
Warning: Explosive substances	20	20	9.105 341
Warning: Toxic substances	20	20	9.105 342
Warning: Corrosive substances	20	20	9.105 343
Warning: Radioactive and ionizing substances	20	20	9.105 344
Warning: Laser beams	20	20	9.105 345
Warning: Oxidizing agents	20	20	9.105 346
Warning: Explosive atmosphere	20	20	9.105 348
Warning: Biodangers	20	20	9.105 347