

4 Drum Sumps, polyethylene

- high resistance to acids, alkalis, oil and other aggressive hazardous materials
- for the storage of pollutants
- available with galvanized or polyethylene grid
- capacity with grid: 2 drums of 60 litres
- capacity without grid: 1 drum of 60 litres
- sump for direct placing onto the ground or on a pallet
- with feet: for safe transport with fork lift truck
- with castors: for safe transport without additional help

asecos GmbH

Description	Grid	Dimensions (W x D x H) mm	Weight kg	PK	Cat. No.
without feet	without	714 x 525 x 235	5	1	9.840 630 1
without feet	galvanized	714 x 525 x 235	12	1	9.840 631
without feet	PE	714 x 525 x 235	13	1	9.840 634
with feet	without	714 x 525 x 330	10	1	9.840 635
with feet	galvanized	714 x 525 x 330	16	1	9.840 633 2
with feet	PE	714 x 525 x 330	17	1	9.840 632
with castors	without	725 x 525 x 380*	15	1	9.840 638
with castors	galvanized	725 x 525 x 380*	20	1	9.840 637 3
with castors	PE	725 x 525 x 380*	22	1	9.840 636

* height without handle



5 Laboratory Trolley clax Mobil comfort

Stable, multi-use trolley with parking brake.
Push-button collapsing/unfolding

Load capacity:

Total: 60kg
Upper shelf: 20kg
Lower shelf: 40kg

Dimensions

unfolded (L x W x H): 900 x 580 x 1020mm
collapsed (L x W x H): 710 x 550 x 190mm



Type	PK	Cat. No.
Laboratory Trolley clax Mobil comfort, incl. 1 clax folding box, 46 litres	1	9.148 100
clax folding box single, 46 litres	1	9.148 099

1


1 Laboratory trolleys

Hard-wearing, sturdy, PE-HD construction resists chemicals and routine knocks and scratches in use, ensuring a long working life.

Langkavel
Advantages:

- stable and high load bearing capacity
- long life
- not electrically conductive
- easy to clean
- hygienic
- high density polyethylene
- easy to assemble without screws
- legs and bases can be recycled
- top tray has integral handles

Type	Colour	Description	Width mm	Length mm	Height mm	PK	Cat. No.
WT 26	beige	2 trays	460	610	660	1	9.148 200
WT 34	beige	3 trays	460	610	860	1	9.148 212
WT 42	beige	3 trays	460	610	1060	1	9.148 216
Cabinet	beige	suitable for WT 26, 34 and 42	460	610	660	1	9.148 203

Further colours available on request.

2


2 Bowl trolleys

NEW!
Gorr GmbH & CO.KG
Trolley

- Blue-grey solid rubber (non-marking) tyres
- Welded frame for 85 or 100 litres plastic bowls
- Dimensions frame (L x W x H): approx. 750 x 490 x 720mm
- Height with red bowl approx. 760mm, with white bowl approx. 840mm
- 4 swivel castors with 75mm diam.
- ball bearing mounted
- powder coated

Bowl, red

- PE-LD
- 85 litres
- Overall Dimensions (L x W x H): approx. 780 x 530 x 350mm

Bowl, white

- PE-LD
- Tapering sides
- 100 litres
- Overall Dimensions (L x W x H): approx. 895 x 530 x 375mm

Description	Colour	Width mm	Length mm	Height mm	PK	Cat. No.
Trolley	black-silver	750	490	720	1	9.148 300
Bowl	red	780	630	350	1	9.148 301
Bowl	white	895	530	375	1	9.148 302

3


3 Trolley, Multiboy

Practical trolley for safe transport in the laboratory.

Stable, two tier, steel frame with push-handle and 4 castors (2 of which can be locked). With 600 x 400mm top shelf and transport box on bottom shelf (capacity 23 litres).

Dimensions (WxDxH):	410 x 660 x 820mm
Height above ground, bottom tray:	130mm
Height above ground, top tray:	500mm
Castor Ø:	75mm
Capacity:	top shelf: 25kg base shelf: 50kg

Type	PK	Cat. No.
Multiboy	1	9.733 670

Trolley, Multiboy, trays

Waterproof plastic tray bottom for the Unilab "Multiboy" trolley. Both sides coated in 4mm (400x600mm) thick melamine resin.

Type	PK	Cat. No.
Trolley, Multiboy, trays	1	9.733 672